

CREATIVE & RESPONSIBLE CHEMISTRY

PRODUCT CATALOGUE

March 2025



INNOVATION
at the heart of
polymers chemistry

THE CHEMISTRY

sustainable and ethical

COP manufactures, synthesizes and formulates a range of materials designed for polymer transformation in several lines of applications, such as orthopedic.



A PEOPLE FOCUSED

Company

French and independent since 1987,
COP develops innovative tailor-made solutions
and can produce in small quantities
to satisfy every specific need.

Composed of chemical engineers and qualified technicians,
the team provides both advice and service.

Certified ISO 9001 v.2015, COP pursues its commitment to
manufacturing to the highest standards..

INNOVATION as a commitment

COP is involved in an intensive research and development
approach and makes available its laboratory for project development.

One goal : to replace toxic components by less harmful equivalent, following
on our first innovations : RTV skin contact silicone range, water expanded
polyurethane foam range, biocomposite : bio-based epoxy resin BPA free and
flax fibers.

Cutting edge equipment : rheometer, tensile testing machine, viscosimeter,
DMA, rheotech, thermoplastic and silicone 3D printer.

Six scientific publications since 2012 : Chemistry of Polymers...

COP is investing in tomorrow's manufacturing process to reduce production
time and costs and meet customization needs. To this end, COP is developing a
range of printable silicones to provide flexible solutions for 3D printing.

“ VISIONARY and precursor

**Because we place people at the heart of our
activities, improving the workplace environment
of our collaborators and partners is our priority.**

This philosophy is also reflected by the conception and
the manufacturing of alternative polymers, which are
bio-based and less dangerous for health.

David Denis, CEO
Since 2010



A large range OF POLYMERS

As a leading expert in chemistry, COP offers its own polymer ranges and customized solutions.

Liquid RESINS

Polymerization yields to materials which characteristics (hardness, flexibility, transparency, elasticity...) make them best suited for various casting techniques.

Depending on their utilization, COP polymers provide : comfort, lightness or mechanical performance.



Silicones

gels, elastomers

Acrylics

composites


Epoxy

biocomposite
(BPA free resin)

Polyurethanes

elastomers, foams,
elastic coating

3



ORTHOPEDICS an historical know-how

For 35 years, COP has been manufacturing and developing new materials for the orthopedics professionals to realize customized equipments.

Since 2016, a training center has been created for complex implementation techniques.

ONE PRIORITY :
Improve comfort for people with disabilities.

Complementary PRODUCTS

COP also supplies all accessories needed for polymers processing operations.

- Molding products and accessories,
- Safety equipment,
- **Natural and standard composite reinforcements.**



Towards TOMORROW'S CHEMISTRY

With its innovative approach as testified by its integrated laboratory, COP is committed to design and produce alternative polymer resins, biosourced and most importantly less harmful for the health of its employees and partners.

A member of the Polyméris competitiveness cluster, it also belongs to a network of academic and private experts working for innovation in the chemistry sector in France.

MATERIALS

Epoxy

COMPOSITE RESINS

- SAFEPOXY@..... 08

RESINS FOR LAMINATION

- ORTHOPOXY@10

Acrylics

RESINS FOR LAMINATION

- COPACRYL.....12

Polyurethanes

HARD FOAMS

- FORMOUSSE14
- ORTHOLEGERE.....15

HARD ELASTOMER

- RESIDUR.....15

SUPPLE FOAMS

- SIPMOUSSE ESTHETIC17
- SIPMOUSSE SUPPLE..... 16

FOAMS IN SHEET

- SIPMOUSSE FINITION.....17

COATINGS

- ERGOPEAU 18
- ERGOFLEX..... 18

Silicones

ELASTOMERS

- COPSIL 3D@.....21
- COPSIL 3D@ ADD-GEL.....22
- COPSIL25
- COPSIL SOCKET27
- COPSIL 6528

GELS

- COPSIL DUPLICATOR.....24
- GEL 00, 25, 30, 40.....24

PASTES

- COPSIL HTV29
- SIPORTHO30

FINISHING COAT

- SILISKIN31

SHEETS

- SILICONE GEL SHEET32

REINFORCEMENTS

ORTHOFLAX REINFORCEMENTS

- Braids.....34
- Non-woven34
- 0/90°34

ECO-BLACK BASALT REINFORCEMENTS

- Braids.....35
- Non-woven35
- 0/90°35

HYBRIDS FLAX/BASALT REINFORCEMENTS

- Braids.....35
- Tape35

CARBON REINFORCEMENTS

- Braids.....36
- Woven mat36

JERSEY & GLASS STOCKINETTE

- Perlon@ stockinette.....36
- Nylglass, 100% glass stockinette 37

ACCESSORIES

LAMINATION

- PVA bags.....38
- Non-woven felt.....38
- Lycra tubes.....38

SILICONE LINER

- Lycra fabric39
- Finishing textiles.....39
- Matrix40
- Distal attachment.....40

CARTRIDGE

- Cartridges, guns, mixing tubes....41

COATING

- Spray gun41

SILICONE CASTING

- RTV silicone casting machine.....42

TECHNICAL PRODUCTS

COLOURING AGENTS

- Colouring agent for PU, Acrylics 44
- Colouring agent for RTV Silicones 44
- Colouring agent for HTV Silicones 44

RELEASING AGENTS

- Releasing agent spray45
- Releasing agent PU grease.....45
- Wax releasing agents46
- Vaseline47
- Releasing agent Z400-V2.....47
- Talc codex.....47

INSULATION AGENTS

- Latex.....48
- Liquid plaster sealer.....48
- Stretch PE film48

GLUES, ADHESIVES & PRIMARY

- Blackglue 0149
- Adhesive tapes49
- Contact glue spray49
- Silicone glue.....50
- Primary grip PM82.....50

SOFTENER

- Gelling agent.....51
- Softening agent51

SOLVENTS


- Acetone52
- Isopropyl alcohol.....52
- Solvent S1 silicone cleaning.....52
- Solvent S3 silicone thinner52

MEASURING & MIXING

- Scales51
- Mixers51
- Pot & buckets51
- Jerrican Tap54

SAFETY & PROTECTION

- Safety suits56
- Gloves56
- Glasses57
- Masks57



TRAINING

Offer

There are a lot of materials in the orthopedic field. Technical knowledge and implementation of polymers are essential in order to obtain performant and reliable devices, meeting the needs of each patient.

That is why **COP Chimie** provides training sessions for orthopedic professionals to discover new materials like **COPSIL 3D®** printable silicones , our flax fibers **ORTHOFLAX®** or in order to perfect your knowledge about a manufacturing technic such as **custom made silicone liner**.

FEW EXAMPLES PROPOSED :

- Silicone 3D printing training with F3DF (module available outside orthopedics for professionals in the industrial, prototyping and luxury sectors)
- Realization of custom made liner with **COPSIL** silicones
- Realization of socket with the new resins **ORTHOPOXY®** and **ORTHOFLAX®** and **ECO-BLACK** reinforcements;
- Realization of tibialis esthetics with **PE SHELLS**
- Etc...



This logo indicates that a training session is available in reference to this material.

► COMPOSITE RESINS

SAFEPOXY®

Epoxy resins for stratification

Workshop sheet

SAFEPOXY® CONTACT

| | |
|---------------|--------------------------------|
| Mixing ratio | 100 / 40 |
| Mix viscosity | 1 100 mPa.s |
| Gel time* | Fast : 15 min Slow : 30 min |

*(measurements made on a 70g mixture at 20°C)

COP makes the DIFFERENCE

SAFEPOXY® resins have lower toxicity compared to market standards.

They are formulated without bisphenol A, an endocrine disruptor identified as SVHC (Substance of Very High Concern), able of interfering with our hormones and producing adverse effects even at very low doses.

Beyond being bisphenol-free and having similar mechanical properties to the marketed epoxy resins, SAFEPOXY® resins are partially bio based. COP succeeded in substituting bisphenol for molecules derived from biomass. The renewable carbon source contained in SAFEPOXY® resins comes from the fermentation of sugars and does not represent any health hazard (INSERM 2016 study).

SAFEPOXY® is our range of epoxy resins dedicated to industrial applications. It is formulated without bisphenol A and from renewable raw materials.

It is aimed at composite processors who want to limit their exposure to hazardous substances and their carbon footprint while manufacturing high-performance composites.

The chemical network is formed and the products are used in the same way as a standard system by mixing the resin and the hardener, applying the resin, crosslinking at room temperature and post-curing to reach the maximum hardness and Tg.

PERFORMANCES

SAFEPOXY® offers the performances of epoxy and combines with both standard fibers (carbon, glass ...) and natural basalt (ECO-BLACK) or linen (ORTHOFLAX®) reinforcements that we offer.

- Transparent products
- High wetting power
- High mechanical tensile and flexural strengths
- Sandable after 16 hours at room temperature
- UV resistant


The different SAFEPOXY® systems reach a Tg close to 80 ° C after post-curing.

We offer a full range of resins to meet different industries' needs.

For laminating

- SAFEPOXY® Contact: contact laminating systems

Laminating systems are available in Slow and Fast reactivities depending on the size of the pieces to be laminated.



► COMPOSITE RESINS

SAFEPOXY®

Epoxy resins for repair and finishing



Workshop sheet

SAFEPOXY® SEALANT

| | |
|--------------|------------|
| Mixing ratio | 100 / 40 |
| Gel time* | 15 min |
| Hardness | 85 Shore D |

SAFEPOXY® GAP FILLER

| | |
|-----------|------------|
| Mix ratio | 100 / 40 |
| Gel time* | 15 min |
| Hardness | 85 Shore D |

*(measurements made on a 70g mixture at 20°C)

For repair

- **SAFEPOXY® Sealant:** consistent paste for repairing or shape refilling

For surface finishing

- **SAFEPOXY® Gap filler:** smoothing and filling coating

RESPECT FOR THE ENVIRONNEMENT

SAFEPOXY® resins are formulated with more than 30% bio-sourced material.

We offer small pack sizes to fit with our customers productions.



| NAMES | | PACKAGING | REFERENCES |
|---------------------|----------------|---------------|------------|
| SAFEPOXY CONTACT | RESIN | 1 kg | SPC R01 |
| | | 2,5 kg | SPC R02 |
| | | 5 kg | SPC R03 |
| | HARDENER SLOW | 400 g | SPC DS01 |
| | | 1 kg | SPC DS02 |
| | | 2 x 1 kg | SPC DS03 |
| | HARDENER FAST | 400 g | SPC DF01 |
| | | 1 kg | SPC DF02 |
| | | 2 x 1 kg | SPC DF03 |
| SAFEPOXY GAP FILLER | RESIN | 1 kg | SPE R01 |
| | | 2,5 kg | SPE R02 |
| | | 5 kg (pot) | SPE R03 |
| | HARDENER | 400 g | SPE D01 |
| | | 1 kg | SPE D02 |
| | | 2 x 1 kg | SPE D03 |
| SAFEPOXY MASTIC | RESIN + FILLER | 1 kg + filler | SPM RC01 |

ORTHOPOXY®

Fast, Slow and Clear



Workshop sheet

ORTHOPOXY® FAST

| | |
|------------------------|--------------|
| Mix viscosity at 20°C | 450 mPa.s |
| Resin / Hardener | 100 / 40 |
| Gel time* | 30 min |
| Gel time in thin layer | 1 h 30 |
| Post curing at 100°C | 1 h |
| Hardness | ± 85 Shore D |

ORTHOPOXY® SLOW

| | |
|------------------------|--------------|
| Mix viscosity at 20°C | 450 mPa.s |
| Resin / Hardener | 100 / 40 |
| Gel time* | 2 h |
| Gel time in thin layer | 5 h 30 |
| Post curing at 100°C | 1 h |
| Hardness | ± 85 Shore D |

ORTHOPOXY® CLEAR

| | |
|------------------------|--------------|
| Mix viscosity at 20°C | 500 mPa.s |
| Resin / Hardener | 100 / 40 |
| Gel time* | 35 min |
| Gel time in thin layer | 3 h 30 |
| Post curing at 100°C | 1 h |
| Hardness | ± 85 Shore D |

*(measurements made on a 70g mixture at 20°C)

New generation of bio-based epoxy resins, BPA free, patented and especially developed for orthopedic equipment.

ORTHOPOXY® range is composed of three versions :

- « FAST » for the realization of prosthetic sockets
- « SLOW » for the realization of orthosis,
- « CLEAR » for the faster realization of aspect parts lightweight and having excellent mechanical and thermal properties. We recommend using the CLEAR version for devices in contact with water or for patients subject to excessive sweating.

These resins have been designed to impregnate our flax fibers (ORTHOFLAX®) and Basalt fibers (ECO-BLACK) but also all the other fibers used in the orthopedic field (CARBON, GLASS, PERLON®, NYLGLASS®, etc.).

- **Comfort at work** : Odorless, Non flammable
- **Performance** : More rigid and more resistant, High compatibility with fibers and easy to impregnate
- **Respect of the environment** : *Bio-based resins*
 - 45% for the FAST version
 - 45% for the SLOW version
 - 50% for the CLEAR version

10



SOCKET IN ORTHOPOXY® AND REINFORCEMENT ORTHOFLAX

| NAMES | | PACKAGING | REFERENCES |
|-----------------|----------|-----------|-------------|
| ORTHOPOXY FAST | RESIN | 1 kg | EPOX-F R01 |
| | | 2,5 kg | EPOX-F R02 |
| | | 5 kg | EPOX-F R03 |
| | HARDENER | 400 g | EPOX-F D01 |
| | | 1 kg | EPOX-F D02 |
| | | 2 x 1 kg | EPOX-F D03 |
| ORTHOPOXY SLOW | RESIN | 1 kg | EPOX-S R01 |
| | | 2,5 kg | EPOX-S R02 |
| | | 5 kg | EPOX-S R03 |
| | HARDENER | 400 g | EPOX-S D01 |
| | | 1 kg | EPOX-S D02 |
| | | 2 x 1 kg | EPOX-S D03 |
| ORTHOPOXY CLEAR | RESIN | 1 kg | EPOX-C R01 |
| | | 2,5 kg | EPOX-C R02 |
| | | 5 kg | EPOX-C R03 |
| | HARDENER | 400 g | EPOX-C D01, |
| | | 1 kg | EPOX-C D02 |
| | | 2 x 1 kg | EPOX-C D03 |

LAMINATING RESINS

COPACRYL

Jersey, carbon, supple, glue, glue gel

Hardener : SIPACRYL

Workshop sheet

JERSEY RESIN

| | |
|-----------------------|--------------|
| Viscosity resin | ± 450 mPa.s |
| Resin/Hardener | 100/2 to 3 |
| Exothermic peak time* | 25 min |
| Demolding time | ± 35 min |
| Hardness | ± 85 Shore D |

CARBON RESIN

| | |
|-----------------------|--------------|
| Viscosity resin | ± 250 mPa.s |
| Resin/Hardener | 100/2 to 3 |
| Exothermic peak time* | 25 min |
| Demolding time | ± 35 min |
| Hardness | ± 85 Shore D |

GLUE RESIN

| | |
|-----------------------|--------------|
| Viscosity | ± 300 mPa.s |
| Resin/Hardener | 100/2 to 3 |
| Exothermic peak time* | 8 min |
| Demolding time | 12 to 15 min |
| Hardness | ± 85 Shore D |

GEL RESIN

| | |
|-----------------------|--------------|
| Viscosity | Thixotropic |
| Resin/Hardener | 100/2 to 3 |
| Exothermic peak time* | 6 min |
| Demolding time | 12 to 15 min |
| Hardness | ± 85 Shore D |

SUPPLE RESIN

| | |
|-----------------------------------|------------------|
| To mix with Jersey /Carbon resins | (10 to 30 %) |
| Viscosity | ± 400 mPa.s |
| Resin/Hardener | 100/2 to 3 |
| Exothermic peak time* | 35 min |
| Demolding time | 35 to 55 min |
| Hardness | 65 to 95 Shore A |

* measurements made on a mixture at 20°C



SOCKET IN COPACRYL CARBON RESIN AND CARBON BRAID

COPACRYL very fluid, quick working time and low odor. Available in Jersey, Carbon, Supple, Glue and Glue-Gel versions.

- Good mechanical properties
- Very good skin tolerance
- Thermoformable after polymerization
(within the limit of elasticity of the fabric used)

New Sipacryl hardener equivalent to the old reference in 2 to 3% dosage, non-CMR.

COPACRYL RANGE

| | |
|-----------------------|---|
| COPACRYL RESIN JERSEY | sockets reinforced with jersey braids |
| COPACRYL RESIN CARBON | sockets reinforced with carbon fiber |
| COPACRYL RESIN SUPPLE | flexibilisation of certain devices (hip socket) in combination with Copacryl Carbon or Jersey |
| COPACRYL RESIN GLUE | liquid and gel |

APPLICATIONS

Can be used with our standard fabrics: **PERLON®**, **NYLON®**, **NYLGLASS®**, **STRETCHNYLGLASS®**, **100% GLASS**.

Can be used with our natural reinforcement **ORTHOFLAX®** and **ECO-BLACK**.

Prosthetic & Orthotics : For lamination of all types of sockets and other splinting devices.

| NAMES | PACKAGING | REFERENCES |
|---------------------------|-------------|-------------|
| COPACRYL JERSEY | 900 g | CAC J01 |
| | 4.9 kg | CAC J05 |
| | 25 kg | CAC J25 |
| COPACRYL CARBON | 900 g | CAC F01 |
| | 4.9 kg | CAC F05 |
| | 25 kg | CAC F25 |
| COPACRYL GLUE | 900 g | CAC C01 |
| | 4.9 kg | CAC C05 |
| | 25 kg | CAC C25 |
| COPACRYL GLUE GEL | 750 g box | CAC G750 |
| COPACRYL SUPPLE | 900 g | CAC S01 |
| | 4.9 kg | CAC S05 |
| | 25 kg | CAC S25 |
| SIPACRYL HARDENER NON-CMR | 150g in bag | SYD-II 101 |
| | 150g in pot | SYDO-II 101 |

- POLYURETHANE

> HARD FOAMS

FORMOUSSE 700, 450, 300, 200

Hard foams more or less expanded

HARDENER MD*

Workshop sheet

FORMOUSSE 700

| | |
|------------------------------|--------------|
| Expansion start time at 20°C | 1 min |
| End of expansion | 2 min |
| Removal from mould | 20 in |
| Expansion | ± 1.4 |
| Density | ± 700 g/L |
| Resin/Hardener | 100/100 |
| Hardness | ± 70 Shore D |

FORMOUSSE 450

| | |
|------------------------------|--------------|
| Expansion start time at 20°C | 1 min |
| End of expansion | 2 min |
| Removal from moul | 10 min |
| Expansion | ± 2.2 |
| Density | ± 450 g/L |
| Resin/Hardener | 100/100 |
| Hardness | ± 55 Shore D |

FORMOUSSE 300

| | |
|------------------------------|--------------|
| Expansion start time at 20°C | 1 min |
| End of expansion | 2 min 20 sec |
| Removal from mould | 10 min |
| Expansion | ± 4.3 |
| Density | ± 230 g/L |
| Resin/Hardener | 100/100 |
| Hardness | 30 Shore D |

FORMOUSSE 200

| | |
|------------------------------|--------------|
| Expansion start time at 20°C | 40 sec |
| End of expansion | 2 min 10 sec |
| Removal from mould | 10 min |
| Expansion | ± 5.9 |
| Density | ± 170 g/L |
| Resin/Hardener | 100/100 |
| Hardness | 25 Shore D |

Range of hard foams more or less expanded (1.5 to 6.5 times), two-component.

All the **FORMOUSSE** foams can be screwed, nailed, stapled and can be worked as wood.

APPLICATIONS IN PROSTHETICS AND ORTHOTICS

FORMOUSSE 200 AND 300 perfectly grip onto acrylic sockets. Depending on its density **FORMOUSSE** foams can be used in orthosis or sockets as a reinforcement in highly stressed areas.

APPLICATIONS IN ORTHOPAEDIC FOOTWEAR

FORMOUSSE 700 and **450** can be used to fill plaster casts. Before filling the cast, we advice to apply either a **LATEX LIQUID** type insulator and then a **RELEASING AGENT Z400** on the inside, or directly a releasing agent such as **DEMOULANT CIRE**.

14

| NAMES | PACKAGING | REFERENCES |
|---------------|-----------|------------|
| FORMOUSSE 700 | 2 kg | FRB 702 |
| | 5 kg | FRB 705 |
| | 60 kg | FRB 760 |
| FORMOUSSE 450 | 2 kg | FRB 402 |
| | 5 kg | FRB 405 |
| | 60 kg | FRB 460 |
| FORMOUSSE 300 | 2 kg | FRB 302 |
| | 5 kg | FRB 305 |
| | 60 kg | FRB 360 |
| FORMOUSSE 200 | 2 kg | FRB 202 |
| | 5 kg | FRB 205 |
| | 60 kg | FRB 260 |
| HARDENER MD* | 2 kg | DMD 002 |
| | 5 kg | DMD 005 |
| | 60 kg | DMD 060 |

* As of August 24, 2023, PROPER TRAINING IS REQUIRED PRIOR TO ANY INDUSTRIAL OR PROFESSIONAL USE.



SHAPES IN FORMOUSSE 450

> HARD FOAMS

ORTHOLEGERE 60

Hard foam, highly expanded

HARDENER MD*

Workshop sheet

| | |
|------------------------------|------------|
| Expansion start time at 20°C | 40 s |
| End of expansion | 2 min |
| Removal from mould | 20 min |
| Expansion | 11.1 |
| Density | 90 g/L |
| Resin/Hardener | 100/100 |
| Hardness | 40 Shore A |



POSITIVE CORSET IN ORTHOLEGERE 60

> HARD ELASTOMER

RESIDUR

Hard elastomer

HARDENER MD*

Workshop sheet

| | |
|-------------------------|------------------|
| Mixing time | 20 / 25 s |
| Gel time | 1 min 30 s |
| Demolding time | 15 min |
| Complete hardening time | 4 h |
| Density | 1 030 g/L |
| Resin/Hardener | 100/45 to 100/55 |
| Hardness | 60 to 70 Shore D |



SHAPE EXTENSION IN RESIDUR

Hard foam highly expanded (12,5 times), two-component. In free expansion or slightly compressed, **ORTHOLÉGÈRE** is designed to fill hollow volumes. This foam can tolerate thermoforming up to 200°C.

Compressed, **ORTHOLÉGÈRE 60** can be used for milling blocks of different density (from 65 to 115 g/l) depending the level of compression.

APPLICATIONS IN PROSTHETICS AND ORTHOTICS

ORTHOLÉGÈRE foam can be used to fill seat- corset buttress or to replace plaster for the realisation of positive seat-corset or orthosis for lower limb.

| NAME | PACKAGING | REFERENCES |
|----------------------|-----------|------------|
| ORTHOLEGERE 60 RESIN | 2 kg | ORT R12 |
| | 5 kg | ORT R15 |
| | 60 kg | ORT R160 |
| HARDENER MD* | 2 kg | DMD 002 |
| | 5 kg | DMD 005 |
| | 60 kg | DMD 060 |

* AS OF AUGUST 24, 2023, PROPER TRAINING IS REQUIRED PRIOR TO ANY INDUSTRIAL OR PROFESSIONAL USE.

Two components hard elastomer unfilled which provides a compact material without bubbles, rigid and hard, which can be flexible depending on the ratio mix used (100/45 : slightly flexible; 100/50 : very slightly flexible ; 100/55 : rigid).

RESIDUR can be screwed, nailed, stapled and can be worked as wood. Slight shrinkage after polymerization : 1,2%.

RESIDUR can also be used as a rigid glue on a lot of rigid materials such as wood (dry).

APPLICATIONS IN ORTHOPAEDIC FOOTWEAR

Used to form the extremity of shapes in FORMOUSSE. After complete cooling, **RESIDUR** is resistant to compression during thermoforming.

| NAMES | PACKAGING | REFERENCES |
|----------------------|-----------|------------|
| RESIDUR RESIN (BLUE) | 2 kg | REDB R02 |
| | 5 kg | REDB R05 |
| HARDENER MD* | 2 kg | DMD 002 |
| | 5 kg | DMD 005 |

* AS OF AUGUST 24, 2023, PROPER TRAINING IS REQUIRED PRIOR TO ANY INDUSTRIAL OR PROFESSIONAL USE.



SIPMOUSSE SUPPLE

Elastic, expanded foam

HARDENER MS*

Workshop sheet

SIPMOUSSE 50 LIGHT

| | |
|------------------------------|-------------|
| Expansion start time at 20°C | 30 sec |
| End of expansion | ± 2 min |
| Removal from mold | 1 h |
| Expansion | 18 |
| Density | 5 g/L |
| Resin/Hardener | 100/80 |
| Hardness | 35 Shore 00 |

SIPMOUSSE 80

| | |
|------------------------------|------------|
| Expansion start time at 20°C | 35 sec |
| End of expansion | ± 3 min |
| Removal from mold | 45 mn |
| Expansion | ± 12 |
| Density | 80 g/L |
| Resin/Hardener | 100/50 |
| Hardness | 5 Shore 00 |

SIPMOUSSE 70/30

| | |
|------------------------------|-------------|
| Expansion start time at 20°C | 35 sec |
| End of expansion | ± 2 min 30 |
| Removal from mold | 45 mn |
| Expansion | ± 12 |
| Density | 80 g/L |
| Resin/Hardener | 100/56 |
| Hardness | 20 Shore 00 |

SIPMOUSSE 30/70

| | |
|------------------------------|-------------|
| Expansion start time at 20°C | 30 sec. |
| End of expansion | ± 3 min |
| Removal from mold | 45 min |
| Expansion | ± 12 |
| Density | 80 g/L |
| Resin/Hardener | 100/60 |
| Hardness | 40 Shore 00 |

Supple and nervous foams, 2 components. Several consistencies are available. After polymerization, **SIPMOUSSE SUPPLE** foams can be easily worked on a router to give a nice surface finish. These foams can be colored and recovered with our elastic and resistant coating : **ERGOPEAU** or **ERGOFLEX**.

APPLICATIONS IN PROSTHETICS AND ORTHOTICS

These foams are used primarily to realize injected cushion seats, in cutom seats and cosmetics.

■ SIPMOUSSE SUPPLE

| | |
|--------------------------------------|--------------------------|
| SIPMOUSSE SUPPLE 50 light | <i>the most expanded</i> |
| SIPMOUSSE SUPPLE 80 supple | <i>the most supple</i> |
| SIPMOUSSE SUPPLE 70/30 medium | <i>intermediary</i> |
| SIPMOUSSE SUPPLE 30/70 firme | <i>the firmer</i> |

| NAMES | PACKAGING | REFERENCES |
|-------------------------------------|-----------|----------------|
| SIPMOUSSE SUPPLE 50 LIGHT RESIN | 2 kg | MSH R02 |
| | 5 kg | MSH R05 |
| SIPMOUSSE SUPPLE 80 SUPPLE RESIN | 2 kg | MSP R02 |
| | 5 kg | MSP R05 |
| SIPMOUSSE SUPPLE 70/30 MEDIUM RESIN | 2 kg | MSL R02 |
| | 5 kg | MSL R05 |
| SIPMOUSSE SUPPLE 30/70 FIRME RESIN | 2 kg | MST R02 |
| | 5 kg | MST R05 |
| HARDENER MS* | 2 kg | DMS 002 |
| | 5 kg | DMS 005 |

* AS OF AUGUST 24, 2023, PROPER TRAINING IS REQUIRED PRIOR TO ANY INDUSTRIAL OR PROFESSIONAL USE

> SUPPLE FOAMS

SIPMOUSSE ESTHETIC

Elastic, firm and lightweight foam

HARDENER MS*

**Workshop sheet**

| | | |
|------------------------------|---|---------------|
| Expansion start time at 20°C | > | 1 min |
| End of expansion | > | 3 min 30 sec |
| Removal from mold | > | 1 h |
| Density | > | ± 80 g/L |
| Resin/Hardener | > | 100/80 |
| Hardness | > | ± 50 Shore 00 |

Supple and nervous foam, two-component which offer an homogenous material, consistent which can support high pressure without totally crushing.

This foam can be colored and recovered by our elastic and resistant coating : **ERGOFLEX**.

APPLICATIONS IN PROSTHETICS AND ORTHOTICS

These foams have been especially designed for the realization of tibialis cosmetics with PE shells.

| NAMES | PACKAGING | REFERENCES |
|--------------------------|-----------|----------------|
| SIPMOUSSE ESTHETIC RESIN | 2 kg | MSC R02 |
| | 5 kg | MSC R05 |
| HARDENER MS* | 2 kg | DMS 002 |
| | 5 kg | DMS 005 |

* AS OF AUGUST 24, 2023, PROPER TRAINING IS REQUIRED PRIOR TO ANY INDUSTRIAL OR PROFESSIONAL USE.


> SHEETS

FOAMS IN SHEET

SIPMOUSSE FINISHING FOAM

Soft foam very elastic in sheets. Use for the recovery of injected custom seats or pressure beds. This foam provides very smooth appearance before the application of our **ERGOPEAU** coating.

| NAME | PACKAGING | REFERENCE |
|--------------------|--------------------|----------------|
| SIPMOUSSE FINITION | 202 x 180 x 0.5 cm | MSP P05 |



ERGOPEAU & ERGOFLEX

Elastic and waterproof coating



Workshop sheet

Spray settings

- 1st layer: round spray, important painting flow
- 2nd layer: smaller spray, horizontal moves, reduced painting flow
- 3rd layer (finishing) : large spray, horizontal moves, important painting flow

Gun pressure

- 1st layer : between 5 et 7 bars
- 2nd layer : 3 bars
- 3rd layer : 5 bars

Drying time at 20° C

ERGOPEAU: 24 h / ERGOFLEX: 1 h

Coloring agents

Do not exceed 2% in weight of the painting (ref. CPU)

Adhesion primer

PU supple foam > without primer

Primer required for closed cells foams rimer

- Plastazote foams type> with primer

18

ERGOPEAU and **ERGOFLEX** are ready-to-use mono-component polyurethane coatings which are very fluid, opaque and which can be colored with our PU color pastes.

These coatings offer a resistant but supple skin, very elastic and waterproof.

ERGOPEAU offers a smooth and glossy skin. It does not shrink after polymerization and 24 hours of drying are needed before manipulating. This coating is adapted to cover beds, seat-corset, for a smooth finish. Before painting, we advise you to cover your support with our finishing foam SIPMOUSSE FINITION in sheet (6mm).

ERGOFLEX offers a matt and stretch skin. Only 1 hour of drying is needed before manipulating and it shrinks a little bit after polymerization which permits to smooth the support (supple). This coating has been especially designed for the finishing and the covering of tibialis or femoral cosmetics.


APPLICATIONS

Protection, finishing, decoration and waterproofing of PU elastomers, rigid or supple.

| NAMES | PACKAGING | REFERENCES |
|------------------|-----------|------------|
| ERGOPEAU COATING | 1 L | EGE3 R01 |
| | 5 L | EGE3 R05 |
| | 25 L | EGE3 R25 |
| | 57 L | EGE3 R57 |
| ERGOFLEX COATING | 1 L | EFX R01 |
| | 5 L | EFX R05 |
| | 25 L | EFX R25 |
| | 57 L | EFX R57 |
| ERGOPEAU PRIMER | 1 L | PPE R01 |
| | 5 L | PPE R05 |
| | 25 L | PPE R25 |
| | 61 L | PPE R61 |

This coating must be used with a vapor aspiration system. To choose the right aspiration system, please contact us.





SILICONE FOR 3D PRINTING

The 3D printing is at the heart of many research themes. The possibilities it offers in terms of design, aesthetics, speed of conception, functionalisation and lightening of structures make it an essential tool for many industries, and in particular for external orthopaedics, which is evolving in the field of personalisation.

In collaboration with F3DF, COP Chimie offers you a training on silicone 3D printing.

> ELASTOMERS

COPSIL 3D®

Printed silicone elastomers



Workshop sheet

COPSIL 3D® 0550

Resin/Hardener 1:5
 Hardness 05 Shore A
 Working time in the mixer* 50 min

COPSIL 3D® 1050

Resin/Hardener 1:5
 Hardness 10 Shore A
 Working time in the mixer* 50 min

COPSIL 3D® 2550

Resin/Hardener 1:5
 Hardness 25 Shore A
 Working time in the mixer* 50 min

COPSIL 3D® 4050

Resin/Hardener 1:1
 Hardness 40 Shore A
 Working time in the mixer* 50 min

* measurements made at 20°C

The **COPSIL 3D®** range is intended for the manufacture of flexible parts of complex design or manufactured in small series (prototyping). The printed silicone elastomers have mechanical characteristics at least equivalent to those of parts made by moulding or injection.

COPSIL 3D® silicone elastomers consist of a two-component system (resin and hardener) mixed in 1:5 parts and 1:1 (**COPSIL 3D 4050**). They are cross-linked at room temperature by polyaddition reaction with a platinum-based catalyst. They are certified for skin contact according to ISO 10993-5.

| NAMES | PACKAGING | REFERENCES |
|----------------|---------------|-------------|
| COPSIL 3D 0550 | KIT 6 x 55 g | 3D0550 S055 |
| | KIT 6 x 850 g | 3D0550 C850 |
| COPSIL 3D 1050 | KIT 6 x 55 g | 3D1050 S055 |
| | KIT 6 x 850 g | 3D1050 C850 |
| COPSIL 3D 2550 | KIT 6 x 55 g | 3D2550 S055 |
| | KIT 6 x 850 g | 3D2550 C850 |
| COPSIL 3D 4050 | KIT 2 x 55 g | 3D4050 S055 |
| | KIT 2 x 850 g | 3D4050 C850 |



COPSIL 3D® ADD-GEL

Silicone printing support



Specificity

COPSIL 3D® ADD-GEL is a support gel for RTV-2 COPSIL 3D® silicone elastomers which benefits are:


- Support complex structures,
- Improve final aspect of printed parts by liquid deposition (smoothing),
- Ready-to-use,
- Transparent (useful for supervising your printing in progress),
- Chemically inert with silicone,
- Water rinseable,
- No impact on mechanical properties of printed silicones,
- Neutral pH,
- No hazard pictogram.

COPSIL 3D® ADD-GEL is a support for RTV-2 liquid silicone printing, allowing the realization of complex parts with important bridges overhangs.

The extrusion of liquid silicone takes place within the gel that supports the printed part. It is ready to use, non-toxic and sinkable.

Its specific consistency allows to smooth the layers of silicone. ADD-GEL is inert and does not impact the mechanical properties of the printed material.

| NAME | PACKAGING | REFERENCE |
|-------------------|-------------|-------------|
| COPSIL 3D ADD-GEL | 1.1 kg pot | ADD-GEL 001 |
| | 5 kg bucket | ADD-GEL 005 |



► GELS

COPSIL GEL 00, 25, 30 & 40

Soft and elastic gels, with fast or slow return

Workshop sheet

COPSIL GEL-00

| | |
|------------------------|------------|
| Working time at 20 °C | 1 h 40 min |
| Demolding time at 20°C | 6 h |
| Final hardening time | 24 h |
| Resin/Hardener | 1:1 |
| Hardness | 0 Shore 00 |

COPSIL GEL-25

| | |
|------------------------|-------------|
| Working time at 20 °C | 1 h 30 min |
| Demolding time at 20°C | 3 h |
| Final hardening time | 24 h |
| Resin/Hardener | 1:1 |
| Hardness | 25 Shore 00 |

COPSIL GEL-30

| | |
|------------------------|-------------|
| Working time at 20 °C | 5 min |
| Demolding time at 20°C | 1 h |
| Final hardening time | 24 h |
| Resin/Hardener | 1:1 |
| Hardness | 30 Shore 00 |

COPSIL GEL-40

| | |
|------------------------|-------------|
| Working time at 20 °C | 8 min |
| Demolding time at 20°C | 50 min |
| Final hardening time | 24 h |
| Resin/Hardener | 1:1 |
| Hardness | 40 Shore 00 |

The silicone gels **COPSIL GEL** are two-component systems (resin and hardener) which are very fluid. They can be mixed in equal parts and cure at room temperature by polyaddition reaction with a platinum catalyst which can be accelerated with heat.

The **COPSIL GEL-30 & 40** are dry, translucent, compact, very soft and resistant materials. Hardness available are 30 and 40 Shore 00. Regarding their viscoelasticity they are classified as «fast return» silicone gels.

The **COPSIL GEL-00** is a translucent, compact and very supple material. Its hardness is non measurable (Shore00) because it is too soft. Regarding its viscoelasticity it is classified as «slow return» silicone gel.

The **COPSIL GEL-25** hardness is 25 Shore 00. It has a sticky feel. It can be used in contact with skin as it respect the ISO 10993-5 standard.

APPLICATIONS

COPSIL GEL-30 or **GEL-40** are primarily dedicated to make supple support as heel cushions, patellar rings pads etc. which do not compress under load.

COPSIL GEL-00 & 25 are primarily dedicated to anti-bedsore sheets.

23

| NAMES | PACKAGING | REFERENCES |
|---------------|---------------------------|-------------|
| COPSIL GEL-00 | KIT 1 KG (500 g + 500 g) | GLCS-00 01 |
| | KIT 10 KG (5 kg + 5 kg) | GLCS-00 10 |
| | KIT 50 KG (25 kg + 25 kg) | GLCS-00 50 |
| COPSIL GEL-25 | KIT 1 KG (500 g + 500 g) | GLC-25 01 |
| | KIT 10 KG (5 kg + 5 kg) | GLC-25 10 |
| | CARTRIDGE (2 x 200 ml) | GLC-25 C400 |
| COPSIL GEL-30 | KIT 1 KG (500 g + 500 g) | GES-30 01 |
| | KIT 10 KG (5 kg + 5 kg) | GES-30 10 |
| | CARTRIDGE (2 x 200 ml) | GES-30 C400 |
| COPSIL GEL-40 | KIT 1 KG (500 g + 500 g) | GES-40 01 |
| | KIT 10 KG (5 kg + 5 kg) | GES-40 10 |
| | CARTRIDGE (2 x 200 ml) | GES-40 C400 |

COPSIL DUPLICATOR

Impression silicone



Workshop sheet:

| | |
|------------------------|------------|
| Working time at 20 °C | 3 min |
| Demolding time at 20°C | 10 min |
| Final hardening time | 5 min |
| Resin/Hardener | 1:1 |
| Hardness | 30 Shore A |

The **COPSIL DUPLICATOR** is a two-component silicone system. Mixed in equal parts and cured at room temperature by polyaddition reaction with a platinum catalyst which can be speeded with heat. After polymerization the material is dry, translucent, compact, very resistant and its hardness is 30 Shore A.

COPSIL DUPLICATOR can be used in contact with the skin as it complies with the ISO 10993-5 standard.

APPLICATIONS

The **COPSIL DUPLICATOR** is mainly used to take impression directly on the skin.



| NAME | PACKAGING | REFERENCES |
|-------------------|--------------------------|---------------|
| COPSIL DUPLICATOR | KIT 1 KG (500 g + 500 g) | DUP 01 |
| | KIT 10 KG (5 kg + 5 kg) | DUP 10 |

► ELASTOMERS

COPSIL

Silicone elastomer

**Workshop sheet:**

| | COPSIL 2 | | COPSIL 3 | COPSIL 5 | | COPSIL 12 | |
|------------------------|-------------|-------------|-------------|-------------|--------|--------------|--------------|
| | normal | fast | fast | normal | fast | normal | fast |
| Working time at 20 °C | 50 min | 20 min | 25 min | 1 h | 25 min | 1 h | 30 min |
| Demolding time at 20°C | 5 h | 2 h | 45 min | 3 h | 45 min | 2 h | 1 h 15 |
| Elongation at break | 480 % | | 870 % | 800 % | | 950 % | |
| Mix viscosity | 7 000 mPa.s | 5 000 mPa.s | 6 000 mPa.s | 5 500 mPa.s | | 11 500 mPa.s | 12 500 mPa.s |
| Hardness | 33 Shore 00 | | 40 Shore 00 | 5 Shore A | | 12 Shore A | |

Specificity

> The whole COPSIL range can be used in contact with skin as it respects the ISO 10993-5 standard.

> Density ± 1

> Resin/Hardener 1:1

Polyaddition silicone range, translucent, two-component and skin contact. **COPSIL** range is composed of silicone more or less fluid which have very good mechanical properties.

Sold in liquid ready to mix with a 1/1 ratio. They are also available in cartridges and can be cast with a Silijet machine.

APPLICATIONS

Designed to make silicone liners, reel cushion, patellar rings, etc.

**CUSTOM LINER IN COPSIL**

| NAMES | PACKAGING | REFERENCES |
|------------------|---------------------------|--------------------|
| COPSIL 2 NORMAL | KIT 1 KG (500 g + 500 g) | T-02SN 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-02SN 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-02SN 50 |
| COPSIL 2 FAST | KIT 1 KG (500 g + 500 g) | T-02SR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-02SR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-02SR 50 |
| | CARTRIDGE (2 x 200 ml) | T-02SR C400 |
| COPSIL 3 FAST | KIT 1 KG (500 g + 500 g) | T-03TR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-03TR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-03TR 50 |
| | CARTRIDGE (2 x 200 ml) | T-03TR C400 |
| COPSIL 5 NORMAL | KIT 1 KG (500 g + 500 g) | T-05TN 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-05TN 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-05TN 50 |
| COPSIL 5 FAST | KIT 1 KG (500 g + 500 g) | T-05TR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-05TR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-05TR 50 |
| | CARTRIDGE (2 x 200 ml) | T-05TR C400 |
| COPSIL 12 NORMAL | KIT 1 KG (500 g + 500 g) | T-12TN 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-12TN 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-12TN 50 |
| COPSIL 12 FAST | KIT 1 KG (500 g + 500 g) | T-12TR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-12TR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-12TR 50 |
| | CARTRIDGE (2 x 200 ml) | T-12TR C400 |

COPSIL

Silicone elastomer

**Workshop sheet:**

| | COPSIL 16 | COPSIL 20 | COPSIL 40 | |
|------------------------|-------------|-------------|-------------|--------------|
| | fast | fast | normal | fast |
| Working time at 20 °C | 30 min | 25 min | 20 min | 10 min |
| Demolding time at 20°C | 50 min | 2 h | 12 h | 3 h |
| Elongation at break | 500 % | 900 % | 250 % | 300 % |
| Mix viscosity | 6 000 mPa.s | 6 000 mPa.s | 45 000mPa.s | 55 000 mPa.s |
| Hardness | 16 Shore A | 20 Shore A | 40 Shore A | |

| NAMES | PACKAGING | REFERENCES |
|------------------|---------------------------|--------------------|
| COPSIL 16 FAST | KIT 1 KG (500 g + 500 g) | T-16SR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-16SR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-16SR 50 |
| | CARTRIDGE (2 x 200 ml) | T-16SR C400 |
| COPSIL 20 FAST | KIT 1 KG (500 g + 500 g) | T-20TR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-20TR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-20TR 50 |
| | CARTRIDGE (2 x 200 ml) | T-20TR C400 |
| COPSIL 40 NORMAL | KIT 1 KG (500 g + 500 g) | T-40TN 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-40TN 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-40TN 50 |
| COPSIL 40 FAST | KIT 1 KG (500 g + 500 g) | T-40TR 01 |
| | KIT 10 KG (5 kg + 5 kg) | T-40TR 10 |
| | KIT 50 KG (25 kg + 25 kg) | T-40TR 50 |
| | CARTRIDGE (2 x 200 ml) | T-40TR C400 |

> ELASTOMERS

COPSIL SOCKET

Silicone elastomer

**Workshop sheet:**

| | |
|------------------------|------------|
| Working time at 20 °C | 30 min |
| Demolding time at 20°C | 1 h |
| Density | 1.1 |
| Resin/Hardener | 1:1 |
| Hardness | 38 Shore A |

Polyaddition silicone, fluid, translucent, two-component and easy to demold. **COPSIL SOCKET** has an extremely high viscosity. This product is used to make a soft socket (in orthopedics), molds that reproduce details perfectly, and resistant parts. Sold in liquid ready to mix with a 1/1 ratio.

APPLICATIONS *Prosthetics*

Designed to make soft sockets, molds or for prototyping.




SOFT SOCKET IN COPSIL SOCKET

| NAME | PACKAGING | REFERENCES |
|---------------|---------------------------|-----------------|
| COPSIL SOCKET | KIT 1 KG (500 g + 500 g) | CSS 01 |
| | KIT 10 KG (5 kg + 5 kg) | CSS 10 |
| | KIT 50 KG (25 kg + 25 kg) | CSS 50 |
| | CARTRIDGE (2 x 200 ml) | CSS C400 |

COPSIL 65

Silicone elastomer very firm

**Workshop sheet**

| | |
|---------------------------|------------|
| Working time at 20 °C | 3 min |
| Removal from mold at 20°C | 20 min |
| Final hardening time | 40 min |
| Density | ± 1.1 |
| Resin/Hardener | 1:1 |
| Hardness | 65 Shore A |

Reinforced polyaddition silicone. Two-component system with a very high hardness, about 65 Shore A. Slightly elastic and very tough. Sold in liquid form ready-to-use

APPLICATIONS

Especially designed to copy temporary sockets.

| NAME | PACKAGING | REFERENCES |
|-----------|--------------------------|-------------------|
| COPSIL 65 | KIT 1 KG (500 g + 500 g) | CF-65SR 01 |
| | KIT 10 KG (5kg + 5 kg) | CF-65SR 10 |



28

SOCKET DUPLICATION IN COPSIL 65

> PASTES

COPSIL HTV

ONE COMPONENT

Workshop sheet

| | |
|----------------------|--------------|
| Post curing | 1 h at 110°C |
| Density | ± 1.1 |
| Hardness | |
| COPSIL HTV 35 | ± 35 Shore A |
| COPSIL HTV 55 | ± 55 Shore A |
| COPSIL HTV 70 | ± 70 Shore A |

**COPSIL HTV**

COPSIL HTV are high temperature vulcanization silicone elastomers with high mechanical properties. They crosslink at heat by polyaddition reaction thanks to a platinum salt based catalyst.

COPSIL HTV are available in three hardnesses :

- 35 Shore A
- 55 Shore A
- 70 Shore A

Silicone pieces from COPSIL HTV are implemented thanks to a calendering machine and crosslink within 1 hour at 110°C only. They can be colored thanks to our range of coloring agents in paste designed for HTV silicones.

APPLICATIONS

Manufacturing of prosthetic (socket, partial amputation, etc.) and orthosis.

| NAME | PACKAGING | REFERENCES |
|------------------------------|-----------|------------------|
| ONE-COMPONENT (1K) | | |
| COPSIL HTV 35 MONO COMPONENT | 2 kg | HTV35 M02 |
| | 5 kg | HTV35 M05 |
| COPSIL HTV 55 MONO COMPONENT | 2 kg | HTV55 M02 |
| | 5 kg | HTV55 M05 |
| COPSIL HTV 70 MONO COMPONENT | 2 kg | HTV70 M02 |
| | 5 kg | HTV70 M05 |

> PASTE

SIPORTHO

Silicone in paste



Workshop sheet

| | |
|---------------------------|------------|
| Mixing time | 20 sec |
| Final hardening time 20°C | 4 min |
| Resin/Hardener | 1:1 |
| Hardness | > |
| SIPORTHO 20 | 20 Shore A |
| SIPORTHO 35 | 35 Shore A |
| SIPORTHO 50 | 50 Shore A |


- **SIPORTHO 20** Soft silicone in paste
- **SIPORTHO 35** Firm silicone in paste
- **SIPORTHO 50** Very firm silicone in paste

Polyaddition silicone in paste, hand mixed, two hardness are available. The hardener is flesh color.

APPLICATIONS

Designed to make orthoplasties.

30



| NAME | PACKAGING | REFERENCES |
|-------------|---------------------------|------------------|
| SIPORTHO 20 | KIT 500 G (250 g + 250 g) | SIP20 500 |
| | KIT 1 KG (500 g + 500 g) | SIP20 01 |
| | KIT 50 KG (25 kg + 25 kg) | SIP20 50 |
| SIPORTHO 35 | KIT 500 G (250 g + 250 g) | SIP35 500 |
| | KIT 1 KG (500 g + 500 g) | SIP35 01 |
| | KIT 50 KG (25 kg + 25 kg) | SIP35 50 |
| SIPORTHO 50 | KIT 500 G (250 g + 250 g) | SIP50 500 |
| | KIT 1 KG (500 g + 500 g) | SIP50 01 |
| | KIT 50 KG (25 kg + 25 kg) | SIP50 50 |

ORTHOPLASTIE IN SIPORTHO 50

> FINISHING COAT

SILISKIN

Soft-touch silicone coating

Workshop sheet

| | |
|-------------------------|--------------------------------|
| Resin/Hardener/Powder | 100/100/30 |
| Mixing time | 1 min |
| Pot life | 36 h |
| Coefficient of friction | 0.9 |
| Consumption | 100 to 150 g by m ² |
| Post cure | 30 min at 100°C |

Slippery finishing coat with soft touch effect for COPSIL silicones RTV, HTV and LSR silicone. SILISKIN coating drastically reduces the coefficient of friction of silicone surfaces. It is extremely flexible and does not stiffen its support.

It is applied in a thin layer and cures with cooking.

APPLICATIONS: research for a soft, slippery, feel of cured silicones, limit dirt build-up



| NAMES | PACKAGING | REFERENCES |
|---|---------------------------|----------------|
| STARTER KIT SILISKIN | 100 g + 100 g + 30 g | SKS 001 |
| SILISKIN COATING | KIT 200 G (100 g + 100 g) | VSK 200 |
| | KIT 1 KG (500 g + 500 g) | VSK 01 |
| | KIT 10 KG (5 kg + 5 kg) | VSK 10 |
| SILISKIN POWDER | 30 g | PSK 030 |
| | 150 g | PSK 150 |
| | 1.5 kg | PSK 015 |
| SPRAY GUN TREND HD + 600 ML GRAVITY FEED CUP | per unit | PIS E01 |
| 600 ML GRAVITY FEED CUP | per unit | GOD 004 |
| METAL TURBINE Ø 45 MM | per unit | AGR 001 |



SILICONE SHEET



SILICONE GEL SHEET

Covered with a PU film.

Transparent silicone gel sheet with a hardness Shore 00 (not measurable due to its softness). Sticky on the surface, elastic and resistant, it is covered on both sides with a thin removable (or not) film of supple polyurethane (25 µm).

APPLICATIONS

Designed to be incorporated in various custom seats or orthoses to avoid pressure sores.

32

| NAMES | PACKAGING | REFERENCES |
|---|-----------|-------------------|
| SILICONE GEL SHEET 400 x 400 x 3 mm | à l'unité | PGS 003 |
| SILICONE GEL SHEET 400 x 400 x 6 mm | à l'unité | PGS 006 |
| SILICONE GEL SHEET 400 x 400 x 10 mm | à l'unité | PGS 010 |
| SILICONE GEL SHEET 400 x 400 x 15 mm | à l'unité | PGS 015 |
| SILICONE GEL SHEET 400 x 400 x 3 mm + ADHESIVE FILM (ONE SIDE) | à l'unité | PGS 003_FA |
| SILICONE GEL SHEET 400 x 400 x 6 mm + ADHESIVE FILM (ONE SIDE) | à l'unité | PGS 006_FA |
| SILICONE GEL SHEET 400 x 400 x 10 mm + ADHESIVE FILM (ONE SIDE) | à l'unité | PGS 010_FA |
| SILICONE GEL SHEET 400 x 400 x 15 mm + ADHESIVE FILM (ONE SIDE) | à l'unité | PGS 015_FA |

- REINFORCEMENT

ORTHOFLAX®

Flax fiber



ORTHOFLAX® reinforcement has been especially developed and patented for orthopedic devices. Flax fiber provides to the composite some properties never reached by the conventional fibers. Can be used with our ORTHOPOXY® resins and also our COPACRYL resins.

APPLICATIONS in Prosthetics and orthotics : Reinforcement for all types of sockets and other splinting devices.



ORTHOFLAX® BRAIDS

Braids knit with flax fibers.

| NAMES | GRAMMAGE | PACKAGING | REFERENCES |
|-----------------------------|----------|-----------|-----------------|
| ORTHOFLAX BRAID - D. 130 mm | 130 g/mL | 5 m roll | B130 005 |
| ORTHOFLAX BRAID - D. 150 mm | 203 g/mL | 5 m roll | B150 005 |

ORTHOFLAX® NON WOVEN

Non woven fabrics made with flax fibers : unidirectional tape, unidirectional fabric (300g/m²), biaxial fabric +/- 45° (350g/m²).


34



| NAMES | WIDTH (MM) | PACKAGING | REFERENCES |
|-----------------------|------------|-----------|----------------|
| ORTHOFLAX FLAT FABRIC | 50 mm | 50 m roll | TAP 050 |
| ORTHOFLAX UD | 1150 mm | 5 m roll | UD 005 |
| | | 10 m roll | UD 010 |
| | | 20 m roll | UD 020 |
| | | 50 m roll | UD 050 |
| ORTHOFLAX BIAXIAL | 1270 mm | 5 m roll | BX 005 |
| | | 10 m roll | BX 010 |
| | | 20 m roll | BX 020 |
| | | 50 m roll | BX 050 |

ORTHOFLAX® WOVEN

Flax fabric, Satin 0-90° (200g/m²) and Serge 2/2 (300g/m²).



| NAMES | WIDTH (MM) | PACKAGING | REFERENCES |
|-------------------------|------------|-----------|----------------|
| ORTHOFLAX 0-90° Natural | 1270 mm | 5 m roll | SA1 005 |
| | | 10 m roll | SA1 010 |
| | | 20 m roll | SA1 020 |
| | | 50 m roll | SA1 050 |
| ORTHOFLAX TWILL 2/2 | 1000 mm | 5 m roll | SE1 005 |
| | | 10 m roll | SE1 010 |
| | | 20 m roll | SE1 020 |
| | | 50 m roll | SE1 050 |




ECO-BLACK

Basalt fibers

Reinforcement range made with basalt fibers. A compromise between carbon and glass fiber. ECO-BLACK reinforcement can be used with our **ORTHOPOXY®** resins and also our **COPACRYL** resins.

APPLICATIONS in **Prosthetics and orthotics** : Reinforcement for all types of sockets and other splinting devices.



ECO-BLACK BRAIDS

Braids knited with basalt fibers.

| NAMES | GRAMMAGE | PACKAGING | REFERENCES |
|-----------------------------|----------|-----------|------------------|
| ECO-BLACK BRAID Ø130mm, 5mL | 244 g/mL | 5 m roll | BB130 005 |
| ECO-BLACK BRAID Ø150mm, 5mL | 488 g/mL | | BB150 005 |

ECO-BLACK TWILL 2/2

Twill 2/2 cloth with basalt fibers.

| NAMES | WIDTH | PACKAGING | REFERENCES |
|---------------------|---------|-----------|----------------|
| ECO-BLACK TWILL 2/2 | 1270 mm | 1 m | SE2 001 |
| ECO-BLACK TWILL 2/2 | | 5 m roll | SE2 005 |

ECO-BLACK NON WOVEN

Unidirectional tape made with basalt fibers.

| NAME | WIDTH | PACKAGING | REFERENCES |
|-----------------------|-------|-----------|----------------|
| ECO-BLACK FLAT FABRIC | 50 mm | 50 m roll | TAP 150 |

HYBRIDS

Flax / Basalt

Reinforcement range made with flax and basalt which permits to combine flax and basalt properties. These reinforcements can be used with our **ORTHOPOXY®** resins and also our **COPACRYL** resins.

APPLICATIONS in **Prosthetics and orthotics** : Reinforcement for all types of socket and other splinting devices.



HYBRIDS BRAIDS

Braids knited with flax and basalt fibers.


| NAMES | GRAMMAGE | PACKAGING | RÉFÉRENCES |
|--------------------------------------|----------|-----------|------------------|
| FLAX/BASALT HYBRID BRAID - D. 130 mm | 185 g/mL | 5 m roll | LB130 005 |
| FLAX/BASALT HYBRID BRAID - D. 150 mm | 370 g/mL | | LB150 005 |
| FLAX/BASALT HYBRID BRAID - D. 225 mm | 461 g/mL | | LB225 005 |

HYBRIDS NON WOVEN

Unidirectional tape made with flax and basalt fibers.

| NAME | DIAMETER | PACKAGING | REFERENCES |
|-------------------------------|----------|-----------|----------------|
| FLAX/BASALT HYBRID TAPE, 50mL | 50 mm | 50 m roll | TAP 250 |

CARBON



CARBON BRAIDS


Carbon braid. These braids are available in 2 diameters (tibial and femoral) and are sold by the kg.

APPLICATIONS in **Prosthetics and orthotics**

Reinforcement of all types of tubular devices.

| NAMES | PACKAGING | REFERENCES |
|-----------------------|---------------------|------------|
| CARBON BRAID D 125 mm | 1 kg (3.7 m at 45°) | GCA 001 |
| CARBON BRAID D 200 mm | 1 kg (3.1 m at 45°) | GCA 002 |

36



WOVEN CARBON MAT

Woven mat at 200 g/m², sold by linear meter.

APPLICATIONS in **Prosthetics and orthotics**

Reinforcement for prosthetic sockets or orthoses.

| NAME | PACKAGING | REFERENCES |
|----------------------|-----------------|------------|
| WOVEN MAT (200 G/M2) | by linear meter | TCA 001 |


JERSEY

Tubular sheaths

Sold by the kilo. 3 diameters available.

APPLICATIONS in **Prosthetics and orthotics**

Reinforcement of all types of tubular devices (rigid sockets).



PERLON® Fibers

| NAMES | PACKAGING | REFERENCES |
|--------------------------|-------------|------------|
| TUBULAR SHEATH - D 10 CM | 1 kg (37 m) | PER 010 |
| TUBULAR SHEATH - D 12 CM | 1 kg (33 m) | PER 012 |
| TUBULAR SHEATH - D 15 CM | 1 kg (27 m) | PER 015 |

GLASS

Tubular sheaths

Sold by the kilo. 7 diameters available.

APPLICATIONS in **Prosthetics and orthotics** : Reinforcement of all types of tubular devices.




NYLGLASS® & STRETCHNYLGLASS® 30% polyamide and 70% glass

STRETCHNYLGLASS® is more elastic than traditional NYLGLASS®. White color.

| NAMES | PACKAGING | REFERENCES |
|---------------------------|-------------|----------------|
| TUBULAR SHEATH - D. 15 CM | 1 kg (22m) | NYL 015 |
| TUBULAR SHEATH - D. 20 CM | 1 kg (18 m) | NYL 020 |


| NAMES | PACKAGING. | REFERENCES. |
|---------------------------|-------------|----------------|
| TUBULAR SHEATH - D. 9 CM | 1 kg (33 m) | SNY 009 |
| TUBULAR SHEATH - D. 12 CM | 1 kg (25 m) | SNY 012 |
| TUBULAR SHEATH - D. 15 CM | 1 kg (21 m) | SNY 015 |



100% GLASS

| NAMES | PACKAGING | REFERENCES |
|---------------------------|------------|----------------|
| TUBULAR SHEATH - D. 15 CM | 1 kg (9 m) | TVE 015 |
| TUBULAR SHEATH - D. 20 CM | 1 kg (7 m) | TVE 020 |

ACCESSORIES FOR LAMINATION



PVA BAGS

Standard PVA bags. 5 sizes are available.

APPLICATIONS

Plaster insulation for sockets or liners lamination.

| NAMES | L X D1 X D2 | PACKAGING | REFERENCES |
|---------------|-----------------|------------|----------------|
| PVA BAG - T.2 | 102 x 5 x 15 cm | Pack of 20 | PVA 015 |
| PVA BAG - T.3 | 102 x 5 x 20 cm | | PVA 020 |
| PVA BAG - T.4 | 102 x 5 x 25 cm | | PVA 025 |
| PVA BAG - T.5 | 102 x 5 x 30 cm | | PVA 030 |
| PVA BAG - T.6 | 102 x 5 x 35 cm | | PVA 035 |

38



NON-WOVEN FELT

Non-woven mat at 200g/m². Available in sheet 2m (1.5m width).

APPLICATIONS in Prosthetics and Orthotics

As an inner layer to resin sockets, this mat gives a very smooth inner surface.

| NAME | PACKAGING | REFERENCE |
|----------------|------------------|----------------|
| NON-WOVEN FELT | sheet 2 m x 1.5m | FEU 001 |




LYCRA TUBES

110 cm length, 20 denier thickness, sold in pack of 50 pieces.

| NAME | PACKAGING | REFERENCE |
|------------------------|------------|----------------|
| TUBULAR LYCRA - 110 CM | pack of 50 | LYT 001 |

ACCESSORIES FOR SILICONE LINER



COVERING MATERIAL

Two way stretch elastic lycra material (polyamide/elastane), in flesh, perfect for covering supple elastomers (depending of material, use the appropriate glue).

APPLICATIONS in Prosthetics and Orthotics

Outer cover for custom made liner.

| NAME | PACKAGING | REFERENCE |
|--------------------|-------------|-----------|
| LYCRA FABRIC FLESH | 2 m x 1,5 m | LYC C01 |

FINISHING TEXTILES

Abrasion-resistant self-smoothing textiles for the manufacturing of custom-made liners. Measuring 40 cm in length, they are available in two formats:

Size S

- Distal circumference: 22 cm
- Proximal circumference: 23 cm

Size L

- Distal circumference: 27 cm
- Proximal circumference: 28 cm

APPLICATIONS in Prosthetics and Orthotics

Custom-made liner manufacturing

| NAMES | PACKAGING | REFERENCES |
|---------------------|-----------|------------|
| FINISHING TEXTILE S | per unit | COV S01 |
| FINISHING TEXTILE L | per unit | COV S02 |



ACCESSORIES FOR SILICONE LINER

MONO-ELASTIC TECHNICAL MATRIX

The matrix support and reinforce the distal attachments. They are mono-elastic distally, elastic proximally, and help counteract the effects of pistoning during the pendulum phase.

APPLICATIONS in Prosthetics and Orthotics

Custom-made liner manufacturing

SHORT (9 CM)

Proximal circumference: 15, 19, 22, 26, 29 and 31 cm.



| NAMES | PACKAGING | REFERENCES |
|---------------------------|-----------|----------------|
| SHORT MATRIX PINK 15 cm | per unit | MAT S01 |
| SHORT MATRIX BLUE 19 cm | per unit | MAT S02 |
| SHORT MATRIX YELLOW 22 cm | per unit | MAT S03 |
| SHORT MATRIX RED 26 cm | per unit | MAT S04 |
| SHORT MATRIX GREY 29 cm | per unit | MAT S05 |
| SHORT MATRIX GREEN 31 cm | per unit | MAT S06 |

LONG (38 CM)

Distal and proximal circumferences:

- 29 and 30 cm
- 32 and 33 cm



| NAMES | PACKAGING | REFERENCES |
|-----------------------------|-----------|----------------|
| LONG MATRIX BLUE 29;30 cm | per unit | MAT L02 |
| LONG MATRIX YELLOW 32;33 cm | per unit | MAT L03 |

DISTAL ATTACHMENTS

For distal attachment liners. Available in 5 diameters. To be drowned into custom made liners.

APPLICATIONS in Prosthetics and Orthotics

For retention of liner to a lock.



| NAMES | PACKAGING | REFERENCES |
|---------------------------------------|-----------|----------------|
| DISTAL ATTACHMENT - SIZE 1 - D. 40 MM | per unit | ATA F01 |
| DISTAL ATTACHMENT - SIZE 2 - D. 50 MM | per unit | ATA F02 |
| DISTAL ATTACHMENT - SIZE 3 - D. 60 MM | per unit | ATA F03 |
| DISTAL ATTACHMENT - SIZE 4 - D. 70 MM | per unit | ATA F04 |
| DISTAL ATTACHMENT - SIZE 5 - D. 80 MM | per unit | ATA F05 |

> CARTRIDGE

CARTRIDGES, GUNS, MIXING TOOLS



COPSIL resins are available in two-component cartridge which allow to inject directly without bubbles.

Empty cartridges and accessories (joints and corks etc.) are available to pack by yourself **COPSIL** resins in cartridge.

Manual and pneumatic two-part gun are also available.

APPLICATIONS

Injection of silicone resins where the mix ratio is 1/1.

| NAMES | PACKAGING | REFERENCES |
|---------------------------------------|------------|------------------|
| CARTRIDGE & ACCESSORIES 2 X 200ML | pack of 5 | CAR 400 |
| 2K DISPENSING GUN PNEUMATIC 2 X 200ML | per unit | PIS CP400 |
| 2K DISPENSING GUN MANUAL 2 X 200ML | per unit | PIS C400 |
| MIXER FOR CARTRIDGE (18 ELEMENTS) | pack of 20 | EMB 818 |
| MIXER FOR CARTRIDGE (24 ELEMENTS) | pack of 20 | EMB 824 |

41

> COATING

SPRAY GUN and accessories



The spray gun has one product entry point at the top (gravity spray gun). The flask proposed is 600mL capacity.

The standard nozzle is 1.6 but can take nozzles up to 3.5. All parts in contact with the liquid are in stainless steel. Cleaning the gun after spraying is easy and flasks are available on request.

APPLICATIONS in Prosthetics and Orthotics

Allows a waterproof and resistant skin on supple foams.

| NAME | PACKAGING | REFERENCES |
|---|-----------|----------------|
| SPRAY GUN TREND HD + 600 ML GRAVITY FEED PAINT CUP | per unit | PIS E01 |
| 600 ML GRAVITY FEED PAINT CUP | per unit | GOD 004 |

RTV SILICON CASTING MACHINE

SILIJET

The machine is equipped with two 2.5 litre product tanks and a manually operated dispensing gun. It is compact and specially designed for small production runs.

The machine uses high accuracy positive displacement piston pumps to provide a constant flow rate for 1:1 such as **COPSIL RTV**.


It is easily connected to the compressed air network..

- Compact and portable
- 5 kg capacity
- More economical and environmentally friendly than cartridges
- Can be sprayed

APPLICATIONS

In Prosthetics and Orthotics : custom-made silicone liner with **COPSIL** silicones.

| NOMS | PACKAGING | RÉFÉRENCES |
|---------------------------|------------|------------|
| SILIJET - CASTING MACHINE | | SIL MC |
| MIXER FOR SILIJET | pack of 20 | EMB 924 |



- TECHNICAL PRODUCTS

► COLORING AGENTS

COLOURING AGENTS

For latex, PU, acrylics and RTV / HTV silicones

Workshop sheet :

- Do not use universal colouring agents.
- Depending on the nature of material use the corresponding colouring agent.
- Maximum concentration: 2% in weight.
- Mix well.

Our range of colouring agents is composed of 3 families :

- Acrylic, PU,
- Silicone,
- Latex.

REMARK

We advice not to put more than 2% in weight in your formulations. At higher level, coloring agent can bring releasing properties and reduce adhesion of coatings such as **ERGOPEAU** or **ERGOFLEX**.

| COLOURING AGENTS FOR PU AND ACRYLIC | | | |
|-------------------------------------|------------|-----------|------------|
| | NAMES | PACKAGING | REFERENCES |
| ● | FLESH | 250 g | CPU C01 |
| ● | DARK FLESH | 250 g | CPU CF1 |
| ● | BLACK | 250 g | CPU N01 |
| ○ | WHITE | 250 g | CPU BL1 |
| ● | RED | 250 g | CPU R01 |
| ● | GREEN | 250 g | CPU V01 |
| ● | YELLOW | 250 g | CPU J01 |
| ● | ORANGE | 250 g | CPU O01 |
| ● | BLUE | 250 g | CPU B01 |

| COLOURING AGENTS FOR RTV SILICONE | | | |
|-----------------------------------|------------|-----------|------------|
| | NAMES | PACKAGING | REFERENCES |
| ● | FLESH | 250 g | CSI C01 |
| | | 500 g | CSI C02 |
| ● | DARK FLESH | 250 g | CSI CF1 |
| | | 500 g | CSI CF2 |
| ● | BLACK | 250 g | CSI N01 |
| | | 500 g | CSI N02 |
| ○ | WHITE | 250 g | CSI BL1 |
| | | 500 g | CSI BL2 |
| ● | RED | 250 g | CSI R01 |
| | | 500 g | CSI R02 |
| ● | GREEN | 250 g | CSI V01 |
| | | 500 g | CSI V02 |
| ● | YELLOW | 250 g | CSI J01 |
| | | 500 g | CSI J02 |
| ● | BLUE | 250 g | CSI B01 |
| | | 500 g | CSI B02 |

| COLORING PASTE FOR HTV SILICONE | | | |
|---------------------------------|------------|-----------|------------|
| | NAMES | PACKAGING | REFERENCES |
| ● | FLESH | 50 g | CSH C50 |
| | | 1 kg | CSH C01 |
| ● | DARK FLESH | 50 g | CSH CF50 |
| | | 1 kg | CSH CF01 |
| ● | BLACK | 50 g | CSH N50 |
| | | 1 kg | CSH N01 |
| ○ | WHITE | 50 g | CSH BL50 |
| | | 1 kg | CSH BL1 |
| ● | RED | 50 g | CSH R50 |
| | | 1 kg | CSH R01 |
| ● | YELLOW | 50 g | CSH J50 |
| | | 1 kg | CSH J01 |
| ● | BLUE | 50 g | CSH B50 |
| | | 1 kg | CSH B01 |

▶ RELEASE AGENTS

RELEASE AGENT FOR HARD PU FOAM

Silicone spray

used with : FORMOUSSE

Workshop sheet

- Shake before use
- Spray in thin layers
- Drying time = 3 min
- Don't touch when the layer is dry.

This silicone spray is a sprayable release agent for hard polyurethane foams and elastomers.

It allows an easy separation from all surfaces except plaster.

Do not forget that nothing adheres on silicone except if you use a primary adhesive.

APPLICATIONS

Allows to remove hard foams from latex. Can be sprayed into PE pots to prevent mixed materials from sticking to the inside of the container. Do not use this releasing agent with **ERGOPEAU**.

| NAME | PACKAGING | REFERENCES |
|------------------------|----------------|------------|
| SILICONE RELEASE AGENT | 500 ml aerosol | ISA 001 |

45

RELEASE AGENT FOR ORTHOLEGERE

PU Grease

Workshop sheet :

- Application by hand using a glove.
- Do not use PU Grease as a separator for flexible foams.
- While the PU Grease is applied don't touch

This paste separator is waxed based, it was developed for separating **ORTHOLEGERE** foam from negative plaster casts.

APPLICATIONS

For separating **ORTHOLEGERE** foams from plaster.

| NAMES | PACKAGING | REFERENCES |
|-----------|-----------|----------------|
| PU GREASE | 5 L | GPU 005 |
| | 5 L | GPU 005 |
| | 50 L | GPU 050 |

> RELEASING AGENT

RELEASING AGENT P 109-V2

for supple PU foams and silicones

Liquid wax based

Liquid wax based releasing agent. For more efficiency, this product can be used filling/emptying inside molds. It leaves an oily film on mold walls which allows an easy demolding.

APPLICATIONS

Can be used with [SIPMOUSSE SUPPLE](#)
Can be used with [COPSIL](#) silicone range.

| NAMES | PACKAGING | REFERENCES |
|---------------------------|-----------|------------|
| WAX RELEASE AGENT P109-V2 | 1 L | DEP-II 001 |
| | 5 L | DEP-II 005 |
| | 56 L | DEP-II 056 |

> RELEASING AGENT

46

RELEASING AGENT P 232-V2

for hard PU foams

Liquid wax based

Workshop sheet

- Mix well the can before use
- Application by brush. Ensure the whole brush is well wetted by the product.
- Drying time at 20° : 40 to 45 mn
- Demoulding will always be easier in making a second layer (wait the first layer is dry).
- Do not touch the film once applied.
- As soon as the foam is hard, still warm, unmold quickly.

Liquid wax based releasing agent. For more efficiency, this product can be used filling/emptying inside molds. It leaves an oily film on mold walls which allows an easy demolding.

APPLICATIONS

Can be used between [FORMOUSSE](#)/plaster.
Can be used between [ORTHOLEGERE](#)/plaster.

| NAMES | PACKAGING | REFERENCES |
|---------------------------|-----------|------------|
| WAX RELEASE AGENT P232-V2 | 1 L | ICI-II 001 |
| | 5 L | ICI-II 005 |
| | 56 L | ICI-II 056 |

REMARKS

Ensure the solvent is completely evaporated before pouring the foam. Dry either in open air or at a temperature not higher than 70°C. Otherwise the foam will have a cratered surface.

It may require several layers if the plaster is thick. In this case ensure the first layer is dry before a second application.

Consider the separator as a layer of paint, it must be thin and smooth to give the desired result.

Do not use this releasing agent on latex, which is attacked by the solvent. In these cases use our Z 400-V2.

> RELEASING AGENT

RELEASING AGENT

for supple PU foams and Acrylics resins

Vaseline

Workshop sheet

- Application by brush in thin layers.
- Do not use Vaseline for rigid foams as FORMOUSSE or ORTHOLEGERE as this will create a soft surface to the foam.
- When the vaseline film is applied, do not touch it anymore

The vaseline is of pharmaceutical quality CODEX and has no risk of irritating the skin. It has the consistency of animal or vegetal fat and is neutral and unalterable in nature, with neither smell or taste. It is a universal release agent for flexible products and some rigid resins and acrylics.

APPLICATIONS

Universal separator for all rigid and flexible products.

| NAME | PACKAGING | REFERENCES |
|----------|-----------|------------|
| VASELINE | 500 ml | VAS 001 |
| | 5 L | VAS 005 |
| | 54 L | VAS 050 |

> RELEASING AGENT

RELEASING AGENT Z 400-V2

for hard PU foams

Liquid silicone based

Workshop sheet

- Application by brush.
- Drying time at 20°C : 40 to 45 min.
- To facilitate separation apply a 2nd layer after the first is dry.
- Don't touch the silicone layer once it is dry
- Remove the piece when the foam is still warm but already hard.

It is made of a solvent which evaporates quickly and in which silicone has been incorporated. It is a universal release agent for all surfaces except plaster. It is very fluid and volatile.

APPLICATIONS

Universal separator for hard foams such **FORMOUSSE** or resins.

| NAME | PACKAGING | REFERENCES |
|---------------------------------|-----------|------------|
| SILICONE RELEASE AGENT Z 400-V2 | 1 L | DEZ-II 001 |
| | 5 L | DEZ-II 005 |
| | 56 L | DEZ-II 056 |

TALC CODEX

Use as filler, sliding agent (liners, PVA bags) or release agent (thermoforming).

| NAME | PACKAGING | REFERENCES |
|------------|-----------|------------|
| TALC CODEX | 500 ml | TAL 001 |
| | 5 L | TAL 002 |
| | 30 L | TAL 030 |

> INSULATORS

INSULATORS

Workshop sheet :

Drying time of the latex skin :

> 20°C : 18 h

48



FOR SUPPLE FOAMS AND PLASTER *liquid latex*

Latex is a milky white product with a strong ammonia smell, which is used as an isolating agent via the thin rubber layer, is left once the water content has evaporated. You can also use the liquid latex by spraying it on fragile supports. The rubber film will then increase their resistance.

APPLICATIONS

In Prosthetics and Orthotics : It is used to manufacture injected seat corset . In this case, latex is sprayed onto the foam which provides resistance and waterproofing. **Podo-orthosis :** Usually used to insulate plaster molds to produce hard foam of feet from FORMOUSSE 450.

| NAME | PACKAGING | REFERENCES |
|------------------------|-----------|----------------|
| ISULATOR LIQUIDE LATEX | 2 L | ILX 002 |
| | 5 L | ILX 005 |
| | 200 L | ILX 200 |

FOR PLASTER

Liquid sealer for dry or wet plaster

LIQUID PLASTER SEALER forms a tight layer that makes it possible to waterproof wet plaster by creating a very thin, smooth, and dry film in just a few minutes. It is ready-to-use and formulated with natural thermoplastic resin.

| NAME | PACKAGING | REFERENCES |
|-----------------------|-----------|----------------|
| LIQUID PLASTER SEALER | 500 ml | IPL 500 |
| | 1 l - 5 l | IPL 001 |
| | 5 l | IPL 005 |



POLYETHYLENE FILM

APPLICATIONS

Prosthetics and Orthotics : To close off plaster casts prior to filled with **ORTHOLEGERE** or **SIPMOUSSE SUPPLE**. **Podo-orthosis :** Separation between patients and plaster casts.

| NAMES | PACKAGING | REFERENCES |
|-----------------|-------------------|----------------|
| STRETCH PE FILM | 150 x 0,10 m roll | FIL 015 |
| STRETCH PE FILM | 300 x 0,45 m roll | FIL 045 |

GLUES & ADHESIVES



BLACKGLUE 01

Non CMR Fast-setting non polyurethane adhesive**

APPLICATIONS Specially formulated for bonding, sealing and repairing a wide range of materials.

| NAMES | PACKAGING | REFERENCES |
|--|---------------------------|--------------------|
| BLACKGLUE 01 NON CMR 50 ml | 1 unit (+ 3 mixing tools) | BG01 C501 |
| BLACKGLUE 01 NON CMR 220 ml | 1 unit (+ 3 mixing tools) | BG01 C221 |
| MIXER FOR 50 ML CARTRIDGES | 3 units | EMBC50 003 |
| | 36 units | EMBC50 036 |
| | 144 units | EMBC50 144 |
| MIXER FOR 220 ML CARTRIDGES | 3 units | EMBC220 003 |
| | 36 units | EMBC220 036 |
| | 144 units | EMBC220 144 |
| 2K DISPENSING GUN MANUAL FOR 50 ML CARTRIDGES | 1 unit | PISC50 001 |
| 2K DISPENSING GUN MANUAL FOR 220 ML CARTRIDGES | 1 unit | PISC220 001 |

49

*According with European Regulation No. 1272/2008, known as CLP for Classification, Labelling and Packaging, as of 13/10/2023.

ADHESIVE PE TAPE

reinforced or non-reinforced

APPLICATIONS Closing molds.

| NAMES | PACKAGING | REFERENCES |
|----------------------------|-------------------|----------------|
| TRANSPARENT NON-REINFORCED | 50 mm x 66 m roll | ADH 001 |
| TRANSPARENT REINFORCED | 50 mm x 66 m roll | ADH A01 |

SUPPLE CONTACT GLUE

APPLICATIONS Versatile adhesive to spray in thin layers.

| NAMES | PACKAGING | REFERENCES |
|-----------------------|----------------|----------------|
| CONTACT GLUE IN SPRAY | 500 ml aerosol | CPA 001 |



► GLUES AND ADHESIVES


GLUES & ADHESIVES

SILICONE GLUE

Mono-component and flexible glue drying with air humidity.

applications Sticking of a lycra fabric on silicone (without using a primary), or silicone on silicone. Can be diluted with our **SOLVENT S3**.

| NAME | PACKAGING | REFERENCES |
|--------------|-------------------------|----------------|
| ACETOXY GLUE | 310 mg cartridge (340g) | COL A02 |



50

ADHESION PRIMER PM82 *for Silicone*

Workshop sheet

| | |
|-------------|----------------------|
| Aspect | Liquid, fluid |
| Density | ± 0.82. |
| Viscosity | 1 mPa.s at 25°C |
| Drying time | mini 15 min, max 6 h |

Grip of the silicone on different supports.

APPLICATIONS The support must be dry, clean and perfectly free of grease. The solvent must be evaporated before pouring the silicone on the support it has to stick with (a white film must be seen).

| NAME | PACKAGING | REFERENCES |
|----------------------|-----------|-------------|
| PRIMARY GRIP PM82 | 250 ml | PM82 |

► SILICON ADDITIVES

GELLING AGENT *for Silicone*

APPLICATIONS It gels RTV silicones, giving it a thixotropic like behavior and allows adjustments when manufacturing custom made liner.

| NAME | PACKAGING | REFERENCES |
|------------------------|-----------|----------------|
| SILICONE GELLING AGENT | 500 ml | AGT 001 |

SILICONE SOFTENING AGENT

APPLICATIONS Silicone oil to reduce silicone resin hardness.

| NAMES | PACKAGING | REFERENCES |
|--------------------------|-----------|-------------------|
| SILICONE SOFTENING AGENT | 2 kg | CNF-50 002 |
| | 5 kg | CNF-50 005 |

► SOLVENTS

SOLVENTS

ACETONE

Solvent and cleaner for the PU and the ERGOPEAU. very flammable, to be used with caution. (Read the safety labels).

APPLICATIONS Stain remover and degreasing agent.

| NAME | PACKAGING | REFERENCES |
|---------|-----------|------------|
| ACETONE | 1 L | ACE 001 |
| | 5 L | ACE 005 |
| | 60 L | ACE 060 |

ISOPROPYL ALCOHOL

Solvent and cleaning agent to be used with caution. Very flammable. (Read the safety labels)

APPLICATIONS Fluid cleaner for PU resins, acrylics and silicones. Stain removal and degreasing agent.

| NAME | PACKAGING | REFERENCES |
|-------------------|-----------|------------|
| ISOPROPYL ALCOHOL | 1 L | ALC 001 |
| | 5 L | ALC 005 |
| | 60 L | ALC 060 |
| | 200 L | ALC 200 |

52

SOLVENT S1 - SILICONE CLEANING

Silicone thinner and cleaner to be used with caution (Read the safety labels).

APPLICATIONS Can be used also as a cleaner for silicones.

| NAME | PACKAGING | REFERENCES |
|------------------------------|-----------|------------|
| SOLVANT S1 SILICONE CLEANING | 1 L | SVS 001 |
| | 5 L | SVS 005 |
| | 60 L | SVS 060 |

SOLVENT S3 - SILICONE THINNER

Silicone thinner to be used with caution (Read the safety labels).

APPLICATIONS Used to dilute silicone glue and very viscous silicone elastomers.

| NAME | PACKAGING | REFERENCES |
|-----------------------------|-----------|------------|
| SOLVANT S3 SILICONE THINNER | 500 ml | SVS3 500 |
| | 1 L | SVS3 001 |
| | 5 L | SVS3 005 |

MEASURING & MIXING




SCALES

- 2 economic models are available depending on the volumes being weighted and the precision of weights required (0,2 kg / 1g), (5kg / 1g).

APPLICATIONS

Precise measuring of chemicals which is essential with PU resins, silicones, acrylics and epoxies.

| NAMES | PACKAGING | REFERENCES |
|---------------------|-----------|----------------|
| SCALE 2000 g to 1 g | per unit | BAL 002 |
| SCALE 5000 g to 1 g | per unit | BAL 003 |



MIXERS

- 3 models in metal for PU foams
- wooden spatula for manual mixing of PU and silicones elastomers, acrylic and epoxies resins

53

APPLICATIONS

Homogeneous mixing of resins and their catalysts.

| NAMES | PACKAGING | REFERENCES |
|-----------------------|------------|----------------|
| METAL TURBINE Ø 45 mm | per unit | AGR 001 |
| METAL TURBINE Ø 65 mm | per unit | AGR 002 |
| METAL TURBINE Ø 90 mm | per unit | AGR 003 |
| WOODEN SPATULA | pack of 50 | SPA B01 |

MIXING CONTAINERS

To obtain a consistent mix we have a series of polyethylene (reusable) pots and buckets in 6 sizes. Their smooth sides improve mixing and removal.



| NAMES | PACKAGING | REFERENCES |
|----------------|------------|----------------|
| POT PP 500 ml | pack of 50 | POT 102 |
| POT PP 1000 ml | | POT 101 |
| BUCKET 2,8 L | pack of 5 | SDM 001 |
| BUCKET 5 L | | SDM 002 |
| BUCKET 17 L | | SDM 003 |
| BUCKET 30 L | | SDM 004 |

MEASURING & MIXING

INJECTION TAP



| NAME | PACKAGING | REFERENCE |
|---------------|-----------|----------------|
| INJECTION TAP | per unit | ROB 006 |

JERRICAN TAP

Small flow tap for 5L or 30L jerrican



| NAMES | PACKAGING | REFERENCES |
|---------------------------|------------|----------------|
| TAP-CAP FOR 30 L JERRICAN | per unit | ROB 004 |
| TAP FOR 5 L CAN | pack of 10 | ROB 007 |

• SAFETY &
PROTECTION

SUIT AND GLOVES



SAFETY SUITS

Disposable lightweight strong work suit. Made from polyethylene (protection category 3, type 4-5-6, to the standard CE 95-0302), it is impermeable to aerosols dust and liquid splashes. It also has a hood and elasticated cuffs and waist and a front zip. Available in 3 sizes : L, XL, XXL.

APPLICATIONS

Protection of people and clothing

| NAMES | PACKAGING | REFERENCES |
|----------------------------------|-----------|------------|
| PROTECTIVE SUIT TIVEK - SIZE L | per unit | COM 001 |
| PROTECTIVE SUIT TIVEK - SIZE XL | per unit | COM 002 |
| PROTECTIVE SUIT TIVEK - SIZE XXL | per unit | COM 003 |

56

NITRILE AND VINYL GLOVES

Nitrile gloves for all kind of product.

⚠ Do not use latex gloves for platinum cured silicones

| NAMES | PACKAGING | REFERENCES |
|-------------------------|-------------|------------|
| NITRILE GLOVES - SIZE M | pack of 100 | GNT 001 |
| NITRILE GLOVES - SIZE L | pack of 100 | GNT L01 |




MASKS AND GLASSES

DUST MASKS

For protection, especially from sanding the rigid and soft elastomers and foams products. Light, hygienic and comfortable, they are held in place with a double elastic band, and are in boxes of 10.

Protection category : EN 149 : 2001 +A1 : 2009 CE0086

| NAMES | PACKAGING | REFERENCES |
|----------------------|------------|----------------|
| DUST PROTECTIVE MASK | pack of 10 | MAS 001 |



VAPOUR MASKS

To protect airways from organic gas or vapours, acid gas or vapours, ammonia and amines. They have a dust proof filter and an active carbon filter, needing no maintenance and disposable once saturated. It is suited for vapour protection from isocyanates from the reaction of our polyurethanes and hardener. These vapours also have suitable warning noticed, and the odour should not be smelt through the mask.

If it is smelly then it is time to change the filter.

The protection class is : FFABEK1P2SL.

| NAMES | PACKAGING | REFERENCES |
|-----------------------|-----------|----------------|
| VAPOR PROTECTIVE MASK | per unit | MAS 002 |



57

PROTECTION GLASSES

Translucid frame glasses for excellent frontal protection, and side screens on the arms. Excellent protection and comfort.

| NAME | PACKAGING | REFERENCES |
|--------------------|-----------|----------------|
| PROTECTIVE GLASSES | per unit | LUN 001 |



INDEX

| | | | |
|------------------------|-------|-------------------------|-------|
| 100% GLASS..... | 38 | NON-WEAVED FELT..... | 36 |
| ADHESIVES..... | 49 | NYLGLASS®..... | 37 |
| ARTGEL COMFORT..... | 32 | ORTHOFLAX®..... | 34 |
| BLAKCGLUE 01..... | 49 | ORTHOLEGERE 60..... | 15 |
| CARTRIDGES..... | 41 | ORTHOPOXY®..... | 10 |
| COATING..... | 19 | PERLON..... | 37 |
| COLORING AGENT..... | 44 | POTS AND BUCKETS..... | 51 |
| COPACRYL..... | 12 | PRIMER..... | 50 |
| COPSIL 65..... | 28 | PROTECTION GLASSES..... | 55 |
| COPSIL DUPLICATOR..... | 25 | PROTECTION MASKS..... | 55 |
| COPSIL GEL..... | 24 | PVA BAGS..... | 39 |
| COPSIL HTV..... | 29 | REINFORCEMENT..... | 34-38 |
| COPSIL SOCKET..... | 27 | RELEASING AGENT..... | 47 |
| COPSIL..... | 26 | RESIDUR..... | 15 |
| COPSIL 3D..... | 22 | SAFETY SUITS..... | 55 |
| COPSIL 3D ADD-GEL..... | 23 | SAFEPOXY®..... | 08 |
| DISTAL ATTACHMENT..... | 40 | SCALES..... | 51 |
| ERGOFLEX..... | 19 | SHEATHS..... | 37 |
| ERGOPEAU..... | 19 | SILESCARE..... | 32 |
| FOAM IN SHEET..... | 18 | SILICONE GLUE..... | 50 |
| FORMOUSSE..... | 14 | SILIJET..... | 42 |
| GELLING AGENT..... | 49 | SILISKIN..... | 31 |
| GLOVES..... | 53 | SIPMOUSSE..... | 17 |
| INSULATOR..... | 48 | SIPMOUSSE ESTHETIC..... | 17 |
| JERRRICAN TAP..... | 51 | SIPMOUSSE FINITION..... | 18 |
| LYCRA MATERIALS..... | 39,40 | SIPORTHO..... | 30 |
| MIXERS..... | 51 | SOFTENING AGENT..... | 51 |
| MIXING JET..... | 39 | SOLVENTS..... | 52 |
| | | SPRAY GUN..... | 41 |
| | | STRETCHNYLGLASS®..... | 38 |
| | | TRAINING..... | 06 |
| | | TWO-PART GUN..... | 41 |
| | | VISCOMOUSSE..... | 18 |



COP

CREATIVE & RESPONSIBLE CHEMISTRY

230 b Route des Bouveries
26190 Saint-Nazaire en Royans
France

Tel 0033 (0) 475 487 720
Fax 0033 (0) 475 450 914

Mail contact@cop-chimie.com

www.cop-chimie.com

