



CREATIVE & RESPONSIBLE CHEMISTRY



MANUFACTURING A CUSTOM LINER

KIT EASYFIT LINER

- ✓ Silicones formulated for liners and adapted to patients' life projects (comfort/proprioception)
- ✓ Includes: 2 x 200 ml silicone cartridges, matrix, distal attachment, glue, finishing textile
- ✓ Technical support for design and manufacturing

Product



230b Route des Bouvieries, 26190 Saint-Nazaire-en-Royans, France



www.cop-chimie.com



contact@cop-chimie.com



+33 (0) 475 487 720

This flyer was designed using the resources of Flaticon.com and printed by Activ' Graphic. - Do not litter

© COP - Version 07-2026

SILICONES IN CARTRIDGES* of 2 x 200 mL

	Product	References	Working time	Demolding time	Hardness
Comfort ↑	COPSIL GEL 25	GLC-25 C400	1 h 30 min	3 h	25 Sh 00
	COPSIL 2	T-02SR C400	20 min	1 h 30	2 Sh A
	COPSIL 3 (white)	T-03TR C400	25 min	45 min	3 Sh A
Proprioception ↑	COPSIL 5	T-05TR C400	20 min	45 min	5 Sh A
	COPSIL 12	T-12TR C400	30 min	1 h 15	12 Sh A
	COPSIL 16	T-16SR C400	30 min	50 min	16 Sh A
	COPSIL 20	T-20TR C400	25 min	2 h	20 Sh A

*2 mixers included

MATRIX

SHORT MATRIX (9 CM)

Reference	MAT S01	MAT S02	MAT S03	MAT S04	MAT S05	MAT S06
Circumference (cm)	15	19	22	26	29	31

LONG MATRIX (38 CM)

Reference	MAT L02	MAT L03
Distal circumference (cm)	29	32
Proximal circumference (cm)	30	33

DISTAL ATTACHMENTS

Reference	ATA F01	ATA F02	ATA F03	ATA F04	ATA F05
Diameter (mm)	Ø 40	Ø 50	Ø 60	Ø 70	Ø 80

FINISHING TEXTILES

Reference	COV S01	COV L01
Length (cm)	50	
Distal circumference (cm)	31.5	33
Proximal circumference (cm)	33	34.5

📍 230b Route des Bouveries, 26190 Saint-Nazaire-en-Royans, France

🌐 www.cop-chimie.com

✉ contact@cop-chimie.com

☎ +33 (0) 475 487 720

This flyer was designed using the resources of Flaticon.com and printed by Activ' Graphic. - Do not litter

© COP - Version 07-2026